

Brain Tumor Diagnosis System Using Vision Transformers

Capstone Project Presentation: Final Requirement for a
Master of Science in Data Science

Andy Achouche, April, 2025.

Objective

- Use a Vision Transformer model to diagnose and classify brain tumors
- Design a web-based application to host the model
- Deploy the application on AWS (Amazon Web Services)
- Secure data in transit between the user and the server

Introduction

- Today's clinical workflows rely heavily on expert radiologists for MRI interpretation, which can be time-consuming and subjective.
- This project introduces an automated system using Vision Transformers to consistently detect and classify brain tumors from MRI scans.
- The solution integrates a state-of-the-art deep-learning model with a user-friendly web interface for seamless image upload and real-time results.
- Fully deployed on AWS with HTTPS/TLS encryption, ensuring secure, scalable access anywhere.

CNN vs. Vision Transformers

- **CNNs:** use convolutional filters for hierarchical feature extraction (local focus)
 - Strengths: efficient parameter sharing, proven performance on **smaller datasets**
 - Limitations: limited global context, potential overfitting on large images
- **Vision Transformers:** split images into patches and apply self-attention across patches (global context)
 - Strengths: capture long-range dependencies, scalable with **large datasets**
 - Limitations: data-hungry, higher computational cost

Data & Preprocessing

- **Dataset Source:** Brain Tumor MRI dataset by Masoud Nickparvar on Kaggle (2021)
- **Composition:** MRI images labeled as glioma, meningioma, pituitary tumor, or no tumor
- **Preprocessing Steps:** skull stripping and bias field correction to remove artifacts and normalize intensities
- **ViT Preparation:** resized images to 224×224 pixels, normalized pixel values, and split into 16×16 patches for transformer input

Model Training & Key Parameters

- **Model Capacity:** dim=512, depth=6, heads=8, mlp_dim=1024
- **Patch & Input:** patch_size=16, image_size=224
- **Regularization:** dropout=0.1, emb_dropout=0.1, weight_decay=1e-5
- **Optimization:** Adam optimizer with learning rate = 1e-4
- **Training Regimen:** batch_size=32, epochs=20, early stopping on validation loss

```
[6]: # Training loop
epochs = 20
loss_values = [] # List to store loss values

for epoch in range(epochs):
    model.train()
    running_loss = 0.0
    for images, labels in train_loader:
        images, labels = images.to(device), labels.to(device)
        optimizer.zero_grad()
        outputs = model(images)
        loss = criterion(outputs, labels)
        loss.backward()
        optimizer.step()
        running_loss += loss.item()


    epoch_loss = running_loss / len(train_loader)
    loss_values.append(epoch_loss) # Store the loss for each epoch

    print(f'Epoch {epoch+1}/{epochs}, Loss: {epoch_loss}')

# Save the trained model
torch.save(model.state_dict(), 'vit_mri_model.pth')
```


Model Performance

- **Accuracy:** 93.4% on the independent test set
- **Confusion Matrix:**
 - Glioma: 266 correct, 33 misclassified (11% error)
 - Meningioma: 268 correct, 38 misclassified (12% error)
 - No Tumor: 404 correct, 1 misclassified (0.2% error)
 - Pituitary: 287 correct, 13 misclassified (4% error)



Web Application Design

- **Front-end:** React single-page interface for image uploads, patient navigation, and results display
- **Back-end:** Flask API routes handling authentication, patient records, test data collection, image inference, and report generation
- **UI Components:** doctor login, patient profile (search/add), medical history, combined tests, imaging instructions, upload, and results pages.

Flask App Architecture

- Built with Flask and SQLite for lightweight, session-based web services
- Modular route design: `/login`, `/patient_profile`, `/medical_history`, `/combined_tests`, `/imaging_studies`, `/upload_mri`, `/predict`, `/generate_report`

Upload MRI Image

Patient: Andy A

Choose File no file selected

Analyze and Classify

Reset

Help

Exit

Doctor Security Code

Enter Security Code:

Submit

Help

Exit

Medical & Family History Assessment

Patient: Andy A

Please answer the following questions:

Medical & Family History

Have you had any prior head trauma? --Select--

Is there any family history of brain tumors? --Select--

Have you experienced unexplained seizures or neurological symptoms? --Select--

Do you have a history of genetic disorders related to tumors? --Select--

Next Step

Help

Exit

Search Existing Patient

Select Patient: --Select a Patient--

Search

Add New Patient

Help

Exit

Authentication & Session Management

- `@app.before_request` guard redirects unauthenticated users to `/login`
- Doctor login via security code; sessions secured with `secret_key` and 10-minute inactivity timeout
- Inactivity timer resets on user events; automatic logout after timeout

Doctor Security Code

Session expires in: 587 seconds

Enter Security Code:

Patient Management Interface

- SQLite `patients.db` stores patient demographics; accessed via `get_db_connection()`
- Add new patient form captures name, age, gender; Search existing patient dropdown populates from the database
- Session stores selected patient for downstream questionnaire and imaging workflows

Add New Patient

Session expires in: 579 seconds

First Name:

Last Name:

Age:

Gender: Male

Submit

HelpExit

MRI Upload & Prediction

- /upload_mri route presents upload form; /predict handles file POST, runs ViT inference
- Preprocessing: resize, normalize; inference with model loaded in-memory
- Outputs: HTML page with uploaded image, histogram of class probabilities, predicted label, and confidence %

Upload MRI Image

Patient: Andy A

Choose File no file selected

Analyze and Classify

Reset

Help


Exit

MRI Classification Result

Session expires in: 581 seconds

Patient: Andy A

Prediction (Histogram):



Predicted Diagnosis: Pituitary Tumor

Confidence: 86.61%

Note: This prediction uses a pre-trained Vision Transformer (ViT) model—a neural network architecture that leverages self-attention mechanisms to capture global image features. Vision Transformers have demonstrated high precision in image analysis and classification.

Upload Another MRI

Start Over (New Patient)

Exit


Generate Medical Report

Help

Report Generation Module

- `/generate_report` compiles patient info, history, test results, MRI findings, and tumor explanation
- Generates a histogram via Matplotlib and embeds it in an HTML report
- Response served with Content-Disposition header for user download as `medical_report.html`

MRI Classification Results



Predicted Diagnosis: Pituitary Tumor
Confidence: 86.61%

Tumor Explanation for Pituitary Tumor

Pituitary tumors occur in the pituitary gland and can affect hormone levels; further evaluation is recommended.

Note: This prediction uses a pre-trained Vision Transformer (ViT) model—a neural network architecture that leverages self-attention mechanisms to capture global image features. Vision Transformers have demonstrated high precision in image analysis and classification.

AWS Deployment & Security

- **Infrastructure:** Dockerized Flask app on EC2 instances behind an Application Load Balancer
- **Scaling:** Auto Scaling group configured for traffic spikes
- **Security:** HTTPS via AWS Certificate Manager; IAM roles least-privilege; security groups restrict access
- **Monitoring:** CloudWatch for logs, metrics, and alerts

Domain & HTTPS Configuration

- **Domain Purchase:** Registered **brain-vit.com** through GoDaddy
- **DNS Management:** Created a hosted zone in AWS Route 53, pointing to Load Balancer
- **SSL/TLS Certificate:** Provisioned via AWS Certificate Manager for brain-vit.com
- **Load Balancer Integration:** Attached certificate to ALB to enforce HTTPS
- **Secure Access:** Application now available at <https://brain-vit.com>, ensuring encrypted data in transit

Future Developments

- **Dataset Expansion:** Incorporate additional imaging modalities and larger, multi-center datasets
- **Edge Deployment:** Optimize model for on-device inference in clinical or mobile settings
- **Continuous Learning:** Implement feedback loops to retrain the model with new real-world cases
- **Data Security Enhancements:** Implement SQLCipher to encrypt the SQLite database at rest for robust application-layer protection
- **Clinical Validation:** Partner with hospitals for pilot studies and regulatory approval

Live Demonstration

Thank you!